

SDK User Manual

Table of Contents

| | |
|---|----|
| 1. Introduction | 3 |
| 2. SDK Overview | 3 |
| 2.1. System Configuration | 3 |
| 2.2. System Requirement | 4 |
| 2.3. API Configuration | 4 |
| 2.4. Data Acquisition Flow chart | 5 |
| 2.5. Acquisition method | 6 |
| 3. API Main function Description | |
| 3.1. ImageCapConnectSensor | 10 |
| 3.2. ImageGetFrameSize | 11 |
| 3.3. ImageCapImageAcquisition | 12 |
| 3.4. ImageCapCommandSend | 14 |
| 3.5. ImageCapStartCapture | 15 |
| 3.6. ImageCapStopCapture | 17 |
| 3.7. ImageCapAcquisitionCancel | 18 |
| 3.8. ImageCapDisConnectSensor | 19 |
| 4. API Sub Function Description | |
| 4.1. ImageCapEnableLog | 20 |
| 4.2. ImageCapGetSAEDSensitivity | 21 |
| 4.3. ImageCapSetSAEDSensitivity | 22 |
| 4.4. ImageCapModeChange | 23 |
| 4.5. ImageCapGetCurrentMode | 24 |
| 4.6. ImageCapGetHostIP | 25 |
| 4.7. ImageCapSearchSensor | 26 |
| 4.8. ImageCapInitSensor | 27 |

| | | |
|-------|---------------------------------|----|
| 4.9. | ImageCapSetWindowTime..... | 28 |
| 4.10. | ImageCapSetAVGModeFrameNum..... | 29 |
| 4.11. | ImageCapGetAVGModeFrameNum..... | 30 |
| 4.12. | ImageCalibration..... | 31 |
| 4.13. | GenerateReferenceFile..... | 33 |

5. Return value Description.....35

6. Reference S/W Description

| | | |
|------|---------------------------|----|
| 6.1. | Overview..... | 36 |
| 6.2. | File Configuration..... | 36 |
| 6.3. | Screen Configuration..... | 36 |
| 6.4. | Sensor Connect | 37 |
| 6.5. | Acquisition..... | 37 |
| 6.6. | Mode Change..... | 38 |
| 6.7. | Window Time | 38 |
| 6.8. | Calibration – Make..... | 38 |

1. Introduction


This document is the SDK User Guide for Bontech Detector. It consists of Network configuration for API, API Usage, and Reference S/W explanations.

2. SDK Overview

Provided APIs are divided into two kinds. The first one is the ImageCaptureDDL which connects with Sensor in order to acquire Images, and the second one is CalibrationDLL for calibration of the acquired Images.

2.1. System Configuration

Detector can be directly connected to PC or indirectly connected through Switch Hub to be used. In connection, Giga-bit Ethernet technology is utilized, and it can transfer up to 1 Gbps using a large amount of Image Data theoretically.



The picture above is an example of System Configuration. Maximum 4 Sensors can be connected to 1 PC.

2.2. System Requirement

- OS : 32bit or 64bit Windows operating systems
- Processor : Dual Core 2.5GHz or faster
- RAM : 1GB or 2GB
- NIC : NIC which support Jumbo Packet

2.3. API Configuration

| | | |
|----------------------------|--|--|
| Header / Library | ImageCapDllEx.h ImageCaptureDLL.lib ImageCaptureDLL.dll GCDLL.dll | <ul style="list-style-type: none"> – Interface with Sensor – Image calibration |
| | CalibrationDllEx.h CalibrationDLL.lib CalibrationDLL.dll | |
| Data File | ImageCapture.ini | –Sensor initialization information file |
| | Ref1.dat | |
| | Ref2.dat | |
| Calibration reference file | [Mode]_Dark.raw | <ul style="list-style-type: none"> – Calibration file by each Mode – 2 types of modes: Trigger mode & AED mode |
| | [Mode]_DefectMap.raw | |
| | [Mode]_Reference00.raw | |

- The function for Detector Interface such as Connection with Sensor, Shooting for Images, and etc. are defined in ImageCapDllEx.h.

When connecting with Detector, you need to set the path of Data file.


- The function for calibrating the X-ray image is defined in CalibrationDllEx.h.

When acquiring the Image, you need to set the path for Calibration reference file.

Calibration reference file can be acquired by using the RawImageViewer program.

- In the Sample S/W, Ref1.dat, Ref2.dat files, and Calibration reference file are not provided.

2.4. Data Acquisition Flow chart



It is a procedure to acquire X-ray images. It is possible to acquire the image if the function is called in order as the flow chart above.

The other functions can be called if necessary.

2.5. Acquisition method

The detector has three methods of image acquisition.

- AED Mode, Soft Trigger Mode, Average Mode, Continuous Mode (Dynamic)

The function to change the mode is ImageCapModeChange.

The image acquisition function is ImageCapImageAcquisition.


The image acquisition function of the Continuous Mode is ImageCapStartCapture.

The following explains when to call ImageCapImageAcquisition and the parameters of the function.

2.6. AED Mode

A method of acquiring images by automatically detecting X-rays.

- Mode change method: ImageCapModeChange (IMAGECAP_AED_MODE);
- Acquisition function: ImageCapImageAcquisition (BRIGHT_MODE);



2.7. Soft Trigger Mode

This is a method of acquiring images from the time of soft trigger occurrence.


- Mode change method:

```
ImageCapModeChange (IMAGECAP_SOFTTRIGGER_MODE);
```

- Acquisition function:

```
ImageCapImageAcquisition (SOFTTRIGGER_INSTANT_MODE);
```

Soft Trigger Mode



- The acquisition function must be called at the same time as the x-ray or after the x-ray.
- The image is acquired at the same time as the function call.

2.8. Average Mode

Average mode is a mode in which the average image is returned after acquiring the set number of frames.

- Mode change method:


```
ImageCapModeChange (IMAGECAP_SOFTTRIGGER_MODE);
```

- Acquisition function:

```
ImageCapImageAcquisition (SOFTTRIGGER_INSTANT_MODE);
```

- Frame number setting method:

```
ImageCapSetAVGModeFrameNum(int nFrameNum);
```



2.9. Continuous Mode (4FPS)

The continuous mode starts by a acquisition command and continuously acquires an image until a stop command is received.

- Mode change method:

```
ImageCapModeChange (IMAGECAP_CONTINUOUS_MODE);
```


- Acquisition function:

```
ImageCapCommandSend (DARK_MODE);
```

```
ImageCapStartCapture (DARK_MODE);
```

```
ImageCapStopCapture ();
```

Continuous Mode



- When to call the Acquisition function : both before and after X-ray irradiation
- Image acquisition starts at the same time as the ImageCapStartCapture function starts
- Acquisition ends when ImageCapStopCapture is called

2.10. Continuous - Binning Mode (8FPS)

The continuous-binning mode starts by a acquisition command and continuously acquires an binning-image until a stop command is received.

- Mode change method:

```
ImageCapModeChange (IMAGECAP_CONTINUOUS_BINNING_MODE);
```


- Acquisition function:

```
ImageCapCommandSend (DARK_MODE);
```

```
ImageCapStartCapture (DARK_MODE);
```

```
ImageCapStopCapture ();
```

Continuous Mode



- When to call the Acquisition function : both before and after X-ray irradiation
- Image acquisition starts at the same time as the ImageCapStartCapture function starts
- Acquisition ends when ImageCapStopCapture is called

3. API Main Function Description

3.1. ImageCapConnectSensor

Header ImageCAPDIIEX.h

Syntax

```
int    ImageCapConnectSensor (
        unsigned int    nHostIP,
        unsigned int    nSensorIP
        LPCTSTR        lpszDataPath
);
```

Parameters

| | |
|---------------------|--------------------------------|
| <i>nHostIP</i> | Host IP address |
| <i>nSensorIP</i> | Sensor IP address |
| <i>lpszDataPath</i> | Sensor Configuration File Path |

Return value

| | |
|-------------|------------|
| IMAGECAP_OK | Success |
| Etc. | Error Code |

Description

By passing the PC IP address and Sensor IP, try to connect to Sensor.

In *lpszDataPath*, ImageCapture.ini, and Ref1.dat, and Ref2.dat files should exist. When the connection is successful, the initialization is performed in *ImagecapConnectsensor*, therefore, image acquisition is possible through ImageCapImageAcquisition() function without calling or configuring for further functions. It does not need to call for ImageCapGetHost() and ImageCapSearchSensor() functions.

It is designed to acquire Image by connecting with Sensor through this function in the Smaple S/W

3.2. ImageCapGetFrameSize

Header ImageCAPDIIEx.h

Syntax

```
int    ImageCapGetFrameSize (  
        unsigned int    nSensorIP,  
        unsigned short* nWidth  
        unsigned short* nHeight  
);
```

Parameters

| | |
|------------------|--------------|
| <i>nSensorIP</i> | Sensor IP |
| <i>nWidth</i> | Frame width |
| <i>nHeight</i> | Frame height |

Return value

| | |
|-------------|------------|
| IMAGECAP_OK | Success |
| Etc. | Error Code |

Description

Obtain the horizontal and vertical size of the frame.

It is 3072x3072 for 4343 sensor, and it is 2500x3052 for 3543 sensor.

3.3. ImageCapImageAcquisition

Header ImageCAPDIIEx.h

Syntax

```
int    ImageCapImageAcquisition (
        unsigned int      nSensorIP,
        unsigned short*    pImage,
        unsigned int      nMode,
        LPCTSTR            lpszRefPath
    );
```

Parameters

| | |
|--------------------|---|
| <i>nSensorIP</i> | Sensor IP |
| <i>pImage</i> | Data buffer to receive Image data |
| <i>nMode</i> | DARK_MODE, BRIGHT_MODE, SOFTTRIGGER_INSTANT_MODE |
| <i>lpszRefPath</i> | Path of Calibration reference file |

Return value

| | |
|-----------------------------|---|
| IMAGECAP_NEW_FRAME | Image data is received. |
| IMAGECAP_RECEIVING | Image data is not received. |
| IMAGECAP_CANCEL | Canceled by a user command |
| IMAGECAP_COMMAND_RECV_ERROR | Data reception error |
| IMAGECAP_RECV_TIMEOUT_ERROR | ACK is not transmitted in a certain time. |
| IMAGECAP_NOTREADY_RO | Not ready to receive Image |
| Etc. | Error Code |

Description

It receives the Image data relevant to *nMode* into *pImage*. *pImage* should be called allocating the size of image data. Since Image size is 16bit, 3072 by 3072, it becomes 18874368byte(3072*3072*(16/8)).

Depending on lpszRefPath configuration, it may obtain calibrated images or non-calibrated images.

| | | |
|-------------|--------------------------|---|
| nMode | DARK_MODE | Acquires a dark image |
| | BRIGHT_MODE | Acquires bright images. When X-ray Expose signal is input, it receives X-ray and then image data. At this time, an infinite waiting for Exposure signal can be cancelled by calling ImageCapAcquisitionCancel() function. Used in AED and Average mode |
| | SOFTTRIGGER_INSTANT_MODE | Used in Soft Trigger Mode. The image is acquired immediately. |
| lpszRefPath | NULL | Acquires un-calibrated images |
| | Calibration path | It acquires calibrated images with reference to file of the inputted path In the Calibration path, the reference file should be prepared in advance. |

Calibration reference files : It is possible to acquire through RawImageViewer program.

- Trigger_Dark.raw, Trigger_DefectMap.raw, Trigger_Referenc00.raw
- AED_Dark.raw, AED_DefectMap.raw, AED_Reference00.raw

3.4. ImageCapCommandSend

Header ImageCAPDIIEx.h

Syntax

```
int ImageCapCommandSend (  
    unsigned int nSensorIP,  
    int nCmd,  
);
```

Parameters

| | |
|------------------|-----------|
| <i>nSensorIP</i> | Sensor IP |
| <i>nCmd</i> | DARK_MODE |

Return value

| | |
|----------------------|----------------------------|
| IMAGECAP_OK | Success |
| IMAGECAP_NOTREADY_RO | Not ready to receive Image |
| Etc. | Socket Error Code |

Description

Transmit command for acquiring images in continuous mode.

When the command is sent, the detector starts to acquire the image.

3.5. ImageCapStartCapture

Header ImageCAPDIIEx.h

Syntax

```
int    ImageCapImageStartCapture (
        unsigned int      nSensorIP,
        unsigned short**  pRefImage,
        unsigned int      nMode,
        LPCTSTR          lpszRefPath
);
```

Parameters

| | |
|--------------------|------------------------------------|
| <i>nSensorIP</i> | Sensor IP |
| <i>pRefImage</i> | Start address of image buffer. |
| <i>nMode</i> | DARK_MODE |
| <i>lpszRefPath</i> | Path of Calibration reference file |

Return value

| | |
|-----------------------------|---|
| IMAGECAP_NEW_FRAME | Image data is received. |
| IMAGECAP_RECEIVING | Image data is not received. |
| IMAGECAP_CANCEL | Canceled by a user command |
| IMAGECAP_COMMAND_RECV_ERROR | Data reception error |
| IMAGECAP_RECV_TIMEOUT_ERROR | ACK is not transmitted in a certain time. |
| IMAGECAP_NOTREADY_RO | Not ready to receive Image |
| Etc. | Error Code |

Description

This function should be called after ImageCapCommandSend.

10 buffers are used internally, and pRefImage receives the starting address of the image being created.

You can acquire a continuous image through the loop statement.

After calling this function, the necessary tasks must be implemented in a thread to maintain the maximum FPS.

It must end the loop statement and call ImageCapStopCapture to end the image acquisition

3.6. ImageCapStopCapture

Header ImageCAPDIIEx.h

Syntax

```
int    ImageCapStopCapture (  
        unsigned int    nSensorIP,  
);
```

Parameters

nSensorIP Sensor IP

Return value

IMAGECAP_OK Success

Description

Image acquisition is stopped in continuous mode.

3.7. ImageCapAcquisitionCancel

Header ImageCAPDIIEx.h

Syntax

```
int    ImageCapAcquisitionCancel (  
        unsigned int    nSensorIP  
);
```

Parameters

| | |
|------------------|-----------|
| <i>nSensorIP</i> | Sensor IP |
|------------------|-----------|

Return value

| | |
|-------------|---------|
| IMAGECAP_OK | Success |
|-------------|---------|

Description

After calling ImageCapImageAcquisition function, cancel the waiting status to receive ready signal. When it is normally canceled, the ImageCapImageAcquisition function will return *IMAGECAP_CANCEL*.

3.8. ImageCapDisconnectSensor

Header ImageCAPDIEx.h

Syntax

```
int    ImageCapDisconnectSensor (  
      unsigned int        nSensorIP  
);
```

Parameters

nSensorIP Sensor IP

Return value

| | |
|--------------------|--|
| IMAGECAP_OK | Success |
| IMAGECAP_STOP_WAIT | Reception thread is not completely finished. |
| Etc. | Error Code |

Description

Close buffer and reception thread allocated in ImageCapConnectSensor() function.

Use only when detector is in single image acquisition mode.

Not used in Dynamic mode.

4. API Sub Function Description

4.1. ImageCapEnableLog

Header ImageCAPDIIEx.h

Syntax

```
int    ImageCapEnableLog (  
      BOOL            nEnable  
);
```

Parameters

| | | |
|----------------|---------|---|
| <i>bEnable</i> | Saving. | TRUE: log is saved FALSE: log is not saved |
|----------------|---------|---|

Return value

| | |
|-------------|---------|
| IMAGECAP_OK | Success |
|-------------|---------|

Description

Set whether to save Log. If you want to save the Log, set nEnable as TRUE. Log files are created according to created dates inside [./log] folder, and as API function is called, the name of function and the return value are printed on log.

4.2. ImageCapGetSAEDSensitivity

Header ImageCAPDIIEx.h

Syntax

```
int    ImageCapGetSAEDSensitivity (
        unsigned int      nSensorIP
        int*              nSens
);
```

Parameters

| | |
|------------------|---|
| <i>nSensorIP</i> | Sensor IP |
| <i>nSens</i> | Receive the current AED sensitivity value |
| | LOW_SENSITIVITY |
| | MID_SENSITIVITY |
| | HIGH_SENSITIVITY |

Return value

| | |
|-------------|------------|
| IMAGECAP_OK | Success |
| Etc. | Error Code |

Description

Receive the currently set AED sensitivity value.

Sensitivity is divided into 3 stages, and the default value is MID_SENSITIVITY.

(LOW_SENSITIVITY, MID_SENSITIVITY, HIGH_SENSITIVITY)

4.3. ImageCapSetSAEDSensitivity

Header ImageCAPDIIEx.h

Syntax

```
int ImageCapSetSAEDSensitivity (
    unsigned int nSensorIP
    int nSens
);
```

Parameters

| | |
|------------------|--|
| <i>nSensorIP</i> | Sensor IP |
| <i>nSens</i> | Sensitivity settings LOW_SENSITIVITY MID_SENSITIVITY HIGH_SENSITIVITY |

Return value

| | |
|-------------|------------|
| IMAGECAP_OK | Success |
| Etc. | Error Code |

Description

Change AED Sensitivity of the sensor. This setting may be one of three values.

LOW_SENSITIVITY, MID_SENSITIVITY, HIGH_SENSITIVITY

4.4. ImageCapModeChange

Header ImageCAPDIEx.h

Syntax

```
int    ImageCapModeChange (  
        unsigned int    nSensorIP,  
        int              nMode  
);
```

Parameters

| | |
|------------------|---|
| <i>nSensorIP</i> | Sensor IP |
| <i>nMode</i> | IMAGECAP_TRIGGER_MODE IMAGECAP_AED_MODE IMAGECAP_SOFTTRIGGER_MODE IMAGECAP_AVG_COUNT_MODE IMAGECAP_AVG_AED_MODE IMAGECAP_CONTINUOUS_MODE |

Return value

| | |
|-----------------|---------------------|
| IMAGECAP_OK | Success |
| IMAGECAP_CANCEL | nMode setting error |
| Etc. | Error Code |

Description

Change the Sensor mode relevant to *nSensorIP*.

4.5. ImageCapGetCurrentMode

Header ImageCAPDIEx.h

Syntax

```
int    ImageCapGetCurrentMode (
        unsigned int    nSensorIP,
        int*            nMode
);
```

Parameters

| | |
|------------------|-------------------------------|
| <i>nSensorIP</i> | Sensor IP |
| <i>nMode</i> | Pointer to receive Mode value |

Return value

| | |
|-------------|------------|
| IMAGECAP_OK | Success |
| Etc. | Error Code |

Description

It receives the current image acquisition method into *nMode*.

nMode – IMAGECAP_TRIGGER_MODE
 IMAGECAP_AED_MODE

If it's not a Success, receive –1 from the nMode.

4.6. ImageCapGetHostIP

Header ImageCAPDIEx.h

Syntax

```
int    ImageCapGetHostIP (
        unsigned int*    pHostIP,
        unsigned int*    pNumberofHost
    );
```

Parameters

| | |
|----------------------|-----------------------------------|
| <i>pHostIP</i> | Pointer of Host IP List |
| <i>pNumberofHost</i> | Pointer of number of Network Card |

Return value

| | |
|------------------------------|---------------------|
| IMAGECAP_OK | Success |
| IMAGECAP_GET_HOST_NAME_ERROR | Get Host Name Error |
| IMAGECAP_NIC_CHECK | No Network cards |
| Etc. | Error Code |

Description

It transmit Network IP List and the number of Network cards of the PC to User Application.

pHostIP has to capture memory as much as the number of Network cards.

4.7. ImageCapSearchSensor

Header ImageCAPDIIEx.h

Syntax

```
int    ImageCapSearchSensor (  
        unsigned int      nHostIP,  
        unsigned int*     pSensorIP,  
        unsigned int*     pNumberOfSensor  
);
```

Parameters

| | |
|------------------------|--|
| <i>nHostIP</i> | IP of Network Card which will check connection |
| <i>pSensorIP</i> | Pointer of Sensor IP connected to nHostIP |
| <i>pNumberOfSensor</i> | Pointer of number of Sensors connected to nHostIP |

Return value

| | |
|-------------|------------|
| IMAGECAP_OK | Success |
| Etc. | Error Code |

Description

It checks whether Network Card for *nHostIP* is connected to Sensor, and when it is connected, it brings IP Address List that is used for communication.

If there are more than one Sensors to Network Card, you can choose which device to use in the User Application.

4.8. ImageCapInitSensor

Header ImageCAPDIEx.h

Syntax

```
int ImageCapInitSensor (  
    unsigned int    nSensorIP,  
    LPCTSTR         lpszDataPath  
);
```

Parameters

| | |
|---------------------|--------------------------------|
| <i>nSensorIP</i> | Sensor IP |
| <i>lpszDataPath</i> | Sensor Configuration File Path |

Return value

| | |
|-------------|------------|
| IMAGECAP_OK | Success |
| Etc. | Error Code |

Description

Refer to the data file in the path of *lpszDataPath*, and initialize the Sensor relevant for *nSensorIP*. The major values are Sensor signal control, M Clock control, Device ID, and etc. *nSensorIP* is the SensorIP found by using ImageCapSearchSensor() function.

Since then, buffer to transmit Image data is internally assigned, and the thread to receive data is practiced. ImageCapDisconnect() function can be summoned and the allocated buffer can be released or the thread can be closed.

4.9. ImageCapSetWindowTime

Header ImageCAPDIEx.h

Syntax

```
int    ImageCapSetWindowTime (  
        unsigned int    nSensorIP,  
        int              nWindowTime  
);
```

Parameters

| | |
|--------------------|---|
| <i>nSensorIP</i> | Sensor IP |
| <i>nWindowTime</i> | Window Time to change. Unit is milliseconds |

Return value

| | |
|-------------|------------|
| IMAGECAP_OK | Success |
| Etc. | Error Code |

Description

This function changes the window time of the detector.

4.10. ImageCapSetAVGModeFrameNum

Header ImageCAPDIEx.h

Syntax

```
int  ImageCapSetAVGModeFrameNum (  
    unsigned int  nSensorIP,  
    int           nFrameNum  
);
```

Parameters

| | |
|------------------|--------------|
| <i>nSensorIP</i> | Sensor IP |
| <i>nFrameNum</i> | frame number |

Return value

| | |
|-------------|------------|
| IMAGECAP_OK | Success |
| Etc. | Error Code |

Description

This function changes the frame number in Average mode.

4.11. ImageCapGetAVGModeFrameNum

Header ImageCAPDIIEx.h

Syntax

```
int  ImageCapGetAVGModeFrameNum (  
    unsigned int    nSensorIP,  
    int*            nFrameNum  
);
```

Parameters

| | |
|------------------|--------------|
| <i>nSensorIP</i> | Sensor IP |
| <i>nFrameNum</i> | frame number |

Return value

| | |
|-------------|------------|
| IMAGECAP_OK | Success |
| Etc. | Error Code |

Description

This function gets the frame number of Average mode

4.12. ImageCalibration

Header

CalibrationDIIEx.h

Syntax

```
int ImageCalibration (  
    unsigned short* pObjectImage,  
    LPCTSTR lpszRefPath,  
    int nWidth,  
    int nHeight,  
    int nMode,  
    CRect rectMargin,  
    BOOL bBackGround  
);
```

Parameters

| | |
|---------------------|---|
| <i>pObjectImage</i> | Object image |
| <i>lpszRefPath</i> | Path of Reference file for Calibration |
| <i>nWidth</i> | Image data width |
| <i>nHeight</i> | Image data height |
| <i>nMode</i> | IMAGECAP_TRIGGER_MODE IMAGECAP_AED_MODE |
| <i>rectMargin</i> | Specify an area to calibrate |
| <i>bBackGround</i> | Set Background processing or not TRUE : Perform Background processing FALSE: Does not perform Background processing |

Return value

| | |
|----|---------|
| 0 | Success |
| -1 | Fail |

Description

Perform Calibration to the X-ray image acquired through Detector.

Necessary reference files are three types by each mode, and the file name should be fixed.

| |
|-------------------------|
| Trigger_Dark.raw |
| Trigger_DefectMap.raw |
| Trigger_Reference00.raw |

If three types of reference files are all existing in the path specified by *lpzRefPath*, perform calibration using these files. The result is written in *pObjectImage*.

You need to call *pObjectImage* after allocating the memory by the size of image size.

If there is no Reference file or if you fail to open the file, -1 will be returned.

Currently this function is not used.

Image correction can be done with ImageCapImageAcquisition().

4.13. GenerateReferenceFile

Header

CalibrationDIIEx.h

Syntax

```
int GenerateReferenceFile (  
    unsigned short*    pDarkImage,  
    unsigned short**   pBrightImage,  
    int                nBrightNum,  
    int                nAcqMode,  
    LPCTSTR            lpszSavePath,  
    int                nWidth,  
    int                nHeight  
);
```

Parameters

| | |
|---------------------|-------------------------|
| <i>pDarkImage</i> | Dark image buffer |
| <i>pBrightImage</i> | Bright image buffer |
| <i>nBrightNum</i> | Number of Bright images |
| <i>nAcqMode</i> | Acquisition mode |
| <i>lpszSavePath</i> | Save path |
| <i>nWidth</i> | Image width |
| <i>nHeight</i> | Image height |

Return value

| | |
|-------|-------------------------------------|
| 0 | No file created |
| Other | Number of generated Reference files |

Description

Generates a reference file for calibration.

Create Gain Map and Defect Map by using Dark image and Bright image.

The file name is different for each mode, and Mode can be checked through ImageCapGetCurrentMode.

The file is saved in SavePath. If this path is passed in during image acquisition, calibrated image acquisition is possible.

5. Return Value Description

Most of the API functions return the IMAGECAP_OK when they succeed.

At the time of Failure, the Error codes are as follows:

| | |
|-------------------------------|--|
| IMAGECAP_SOCKET_LOAD_ERROR | Socket Load Error |
| IMAGECAP_INVALID_SOCKET_ERROR | Socket Formation Error |
| IMAGECAP_BIND_SOCKET_ERROR | Socket Bind Error |
| IMAGECAP_SOCKET_OPTION_ERROR | Socket Option Configuration Error |
| IMAGECAP_COMMAND_SEND_ERROR | Command Transmission Error |
| IMAGECAP_COMMAND_RECV_ERROR | Data Reception Error |
| IMAGECAP_RECV_TIMEOUT_ERROR | ACK is not transmitted in certain time. |
| IMAGECAP_ACK_SEND_ERROR | ACK Transmission Error |
| IMAGECAP_NO_SENSOR_ERROR | No Sensor in Lan area |
| IMAGECAP_REGISTER_TYPE_ERROR | A Wrong Register Type |
| IMAGECAP_GET_HOST_NAME_ERROR | Get Host Name Error |
| IMAGECAP_NIC_CHECK | No Network Cards |
| IMAGECAP_STOP_WAIT | Receiving thread is not completed. |
| IMAGECAP_OK | Success |
| IMAGECAP_CANCEL | Canceled Image Acquisition |
| IMAGECAP_NEW_FRAME | Acquisition of a New Image |
| IMAGECAP_RECEIVING | Failed in Image Acquisition |
| IMAGECAP_NOTREADY_RO | Tried Image Acquisition within 10 sec Image Cycle Time |
| Undefined Error Code | Same as Windows System error code |

6. Reference S/W description

6.1. Overview

It is a simple MFC dialog which can acquire image data and perform calibration using the acquired data through Detector.


6.2. File Configuration

The files necessary for practice consist of 7 files, and the files necessary for Calibration consist of 3 files.

| | | |
|----------------|---------------------|--|
| EXE file / DLL | ImageCAP_sample.exe | EXE file |
| | ImageCaptureDLL.dll | It must be located in the same path as EXE file. |
| | CalibrationDLL.dll | |
| | GCDLL.dll | |
| Data file | ImageCapture.ini | Data file necessary for EXE file. |
| | Ref1.dat | The path should be set when Sensor is connected. |
| | Ref2.dat | |

- Ref1.dat, Ref2.dat : It is possible to execute without these in the sample SW, however, these are required for the actual Sensor test.
- Calibration reference file : Without this file, Non-calibrated images will be acquired. It is not provided in the sample SW

6.3. Screen Configuration



6.4. Sensor Connect

When the connect button is clicked, ImageCapConnectSensor function is called to connect with the detector.

m_strReferencePath is set to ./A_Data.

You should have an ImageCapture.ini file inside the folder.

Then check the FrameSize.

6.5. Acquisition (AED, Soft Trigger, Average Mode)

- Get Dark

This button acquires a dark image and saves it in the A_Cal folder.

The file name is different for each mode.

The file in A_Cal is used to calibrate the X-ray image.

Allocating memory to receive the dark image and using the ImageCapImageAcquisition () function to obtain the image.

If the function returns IMAGECAP_NEW_FRAME, the acquisition was successful.

- X-ray

This button acquires an X-ray image. It is saved as [Image.raw] file in the same path as the executable file.

It works differently depending on the third argument of ImageCapImageAcquisition.

| Acquisition Mode | Argument of ImageCapImageAcquisition |
|------------------|--------------------------------------|
| AED | BRIGHT_MODE |
| Soft Trigger | SOFTTRIGGER_INSTANT_MODE |
| Average | BRIGHT_MODE |

AED mode recognizes X-ray automatically, so if you call Acquisition function, it will be in infinite standby. To cancel this, you need to call ImageCapCancel.

The reference file for calibration was set to [. \ A_Cal] folder.

6.6. Acquisition (Continuous Mode & Continuous Binning Mode)

- Get Dark

The dark image is the same as other mods.

However, 10 dark images were obtained for calibration.

- X-Ray

If you click the x-ray button, Continuous mode starts through ImageCapCommandSend, and Image data receives the address value through ImageCapStartCapture function and copies the memory.

You can stop the acquisition through the cancel button.

6.7. Mode Change

The acquisition mode changes to the selected mode.

6.8. Window Time

The window time changes to the time entered.

6.9. Calibration – Make

Obtain a dark image through the Get Dark button, obtain a bright image through the X-ray button, and pass the buffer to the GerateReferenceFile function to create a calibration file.

For samples s/w, the generated calibration file is stored in A_Cal.

The file name is stored differently depending on the acquisition mode.